

Vital patient data at a glance



Sentec Connectivity enables the transmission of oxygenation and ventilation monitoring - **tcPCO₂, tcPO₂, SpO₂, Pulse Rate (PR), Pulsation Index (PI), Heating Power (HP)** - to the Dräger Integrated Care Manager (ICM).






Sentec Connectivity Dräger Integrated Care Manager (ICM)

- Displays and stores near real-time data from the Sentec Digital Monitor
- The Sentec Digital Monitor is automatically identified when connected to the Integrated Care Manager
- Integrated Care Manager supports and provides a dedicated and ergonomic clinical workflow for clinicians

Care with
Confidence

Sentec Connectivity

Dräger Integrated Care Manager (ICM)

Required Components	System Requirements	Description
Sentec Digital Monitor (SDM) 	Software version V08.00 to V08.03	Provides noninvasive and continuous real-time monitoring of ventilation and oxygenation
Integrated Care Manager (ICM) 	Version 9.2	ICM is a medical device software with Device Class IIa. For further detailed information contact your Dräger representative Displays and stores customizable data from SDM
Standard RS232 cable or USB2Serial converter cable 		Connects the SDM with the PC on which ICM is running

Sentec Connectivity Support Team

Let's connect.

If you need support, please contact your local Sentec representative, send an e-mail to connectivity@sentec.com or consult your Dräger representative.

Please visit www.sentec.com for more information.

Your local distributor:

sentec.

Ringstrasse 39
4106 Therwil
Switzerland
www.sentec.com
 Made in Switzerland

Sentec Connectivity Support Team:
Tel +41 61 726 97 60
connectivity@sentec.com

Dräger

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23558 Lübeck
Deutschland
Tel +49 451 882 0
www.draeger.com

Note: This is a global overview and may contain information on products not released in all markets/countries. Information is subject to change without notice.