

Patient Monitoring Systems (PMS)		
Ascom	Digistat® Smart Central	
Dräger	Infinity® Delta (XL), Infinity® Kappa, Infinity® Omega Serie	
Fukuda Denshi	DS-8400	
GE Healthcare	Carescape B450/B650/B850, Carescape Canvas, Solar 8000 M/i, Dash 5000	
Mindray	BeneView T5/T6/T8/T9, BeneVision N12, N15, N17, N19, N22	
Philips	IntelliVue MP40/50/60/70/80/90, MX400/450/500/550/600/700/750/800/850	
Spacelabs Healthcare	Flexport System Ultra View SL (90442A-03)	
Patient Data Management Systems (PDMS)		
Capsule (Coming Soon)	Medical Device Information Platform (MDIP)	
Dräger	Integrated Care Manager (ICM)	
Fukuda Denshi	Mirrel	
iMDsoft	MetaVision	
Nihon Kohden	Prime Gaia	
Stowood	Visi-Download	
TZM	TZM (via) UMG to GE Healthcare's patient monitors	
Ventilators		
Air Liquide	Vendom 30/40	
Breas	Vivo 45 (LS/L), Nippy 4+	
EOVE	EO-150	
Philips	Trilogy Evo (O2), Trilogy EV300, A30 EFL, A40 Pro/EFL	
Thora Tech	Monsoon 4	
Digital Connectivity to Polygraphic and Polysomnographic Systems (PG/PSG)		
Cidelec	CID102L, CID102LB, CID-LXe, CID-LXe-206d, CID-LX, CID-LX-206d	
Nox Medical	NOX A1, NOX A1s (via C1), NOX C1	
Analog Connectivity to Polygraphic and Polysomnographic Systems (PG/PSG)		
Manufacturer/Brand		PSG Adapter Cable (REF)
Carefusion	SomnoStar Pro	PSG Cable E (103678)
Cidelec	CID102L8-D	PSG Cable P (103686)
Compumedics	E-Series, Grael (4k), Somte	PSG Cable E (103678)
	Siesta	PSG Cable G (103680)
Dr. Fenyves & Gut /Löwenstein Medical	SleepDoc Porti 7/8, MiniScreen PRO/plus	PSG Cable D (103677)
	Sonata	PSG Cable B (103675)
Natus	Comet / Aurora, Embla SDx, Embla Titanium (Isolation Box)	PSG Cable E (103678)
	Comlab, Deltamed Duo 44	PSG Cable H (103681)
	Deltamed Brain Box	PSG Cable J (103682)
	Embletta Gold, Embletta MPR PG	PSG Cable L (103684)
	Embla NDx, PAIB2 interface box	PSG Cable F (103679)
	Embla S7000 / N7000	PSG Cable A (103674)
	Embla Titanium	PSG Cable K (103683)
Mdrive, TX Proxy	PSG Cable N (103685)	
Neurovirtual	Sleepvirtual, BWII (SI1089), BWIII	PSG Cable E (103678)
SOMNOmedics	ADC Board	PSG Cable E (103678) (Note: Connection with the SDM may require additional safety measures to be compliant with the requirements of IEC 60601-1. Refer to SDM Technical Manual)
	SOMNO HD (NGA100)**	PSG Cable Q (103687)
	SOMNO HD eco (NGA100)**	PSG Cable Q (103687)
	SOMNOscreen (plus) (SSC082)	PSG Cable Q (103687)
	SOMNOtouch RESP (TOS120)**	PSG Cable Q (103687)
Stowood Scientific Instruments	Visi-3, Visi-System	PSG Cable E (103678)
Philips Respironics	Alice 3 / 4 / 5 / 6 (square connector)	PSG Cable B (103675)
	Alice 3 / 4 / 5 / 6 (round connector)	PSG Cable F (103679)
Weinmann MCC	SOMNOlab, SOMNOlab2	PSG Cable C (103676)



For a detailed overview of the SDM connectivity options including compatibility, technical information and instructions, refer to sentec.com/transcutaneous-connectivity-overview/

Note: This is a global overview and may contain information on products not released in all markets/countries. Information is subject to change without notice.

Sentec AG
Ringstrasse 39
CH-4106 Therwil
Switzerland
sentec.com

** Product number of accessory cables or optocoupler required for installation in brackets behind the product name.